


PREPARATION

Q-Kit is designed To prepare 2 x 5 ml of Autologous Platelet Rich Plasma (A-PRP). It's recommended for Mesotherapy, fine wrinkles and Alopecia treatments.

- Blood collected from the patient in special tubes (10 ml of blood)
- centrifuged at 3200 RPM for 5 min in order to obtain the PRP The separation of platelets in plasma from red blood cells is allowed by the innovative gel placed into the tubes. The PRP is then ready for its application either alone or mixed with an activating component.



Before



After



Before




After

Applications

IMPROVEMENT IN SKIN QUALITY

- THICKNESS
- COLOR
- TEXTURE
- HYDRATION

B SPECIFIC LOCAL IMPROVEMENT





Autologous Cellular Rejuvenation Platelet Rich Plasma (PRP)

Q-kit PRP

Q-kit PRP is a revolutionary method, based on using ingredients of your own blood, which can regenerate your skin. These ingredients reside within the fraction of your blood named "Platelet Rich Plasma" (PRP). PRP is enriched with platelets stem cells and growth factors. Once injected into your skin these substances will accelerate the normal healing process via new collagen formation at the injected area. This process brings to long lasting tissue rejuvenation by your own cells. PRP has been widely used for many years in plastic surgery, orthopedics, dentistry, soft tissue implants and other types of implants



ADVANTAGES

Simplicity: Only few straightforward steps are required to prepare 4 to 5 ml of Platelet-Rich Plasma (PRP).

Low Blood Volume Required: PCLQ-Kit is designed to prepare Autologous Platelet Rich Plasma from only 10 ml of blood.

Optimum Platelet Recovery:

With PCLQ-Kit you recover around 80% of the platelets from the whole blood. Platelets are concentrated in 1 to 5 ml of autologous plasma (depending on the required platelet concentration).

